

Suture Brand Name Cross Reference Chart

Material	UNIFY®	Ethicon®	Covidien®	B.Braun®
Coated or uncoated braided synthetic mid-term absorption	PGA	Vicryl® Coated	Polysorb® / Dexon® II	Safil®/Novosyn®
Monofilament synthetic mid-term absorption	PGCL	Monocryl®	Biosyn®	Monosyn®
Coated or uncoated braided synthetic fast term absorption	RPGA	Vicryl® Rapide	-	Safil® Quick
Monofilament synthetic long term absorption	PDO	PDS® II	Maxon®	Monoplu
Non Absorbable monofilament polyamide	NYLON	Ethilon®/Ethicrin®	Monosof® /Dermalon®	Dafilon
Non Absorbable monofilament Polypropylene	POLYPROPYLENE	Prolene®	Surgipro®/Surgilene®	Premilene®
Non absorbable braided natural silk	SILK	Perma-Hand®	Sofsilk®	Silkam®
Non absorbable braided polyester (coated)	POLYESTER	Ethibond®	Surgidac®	Flexidene®
Non absorbable multifilament polyamide	CABLE NYLON	Nurolon®	Surgilon®	-